

Novel Ionizable Lipids for mRNA Delivery in Rodents and Non-Human Primates

FUJIFILM

Shinichi Hashimoto¹, Yuki Imaizumi¹, Toshifumi Kimura¹, Daisuke Nakagawa¹, Kohei Yasuda¹, Kohei Shimizu¹, Akira Inomata², Masaki Noro¹, Hirofumi Fukunaga¹, Shigetomo Tsujihata^{1,2}
¹FUJIFILM Corporation ²FUJIFILM Toyama Chemical Co., Ltd.

13th international mRNA health conference
(2025/11/11-11/13, Berlin)

Abstract

[Introduction] Among the constituents of lipid nanoparticles (LNPs), ionizable lipids play a critical role in determining the efficacy and safety of mRNA-LNP therapeutics. We offer comprehensive services for DDS therapeutics, from formulation design to commercial production. We also license our proprietary ionizable lipids. To date, we have identified **FL-1245T** for hepatic delivery that when formulated into LNPs, exhibit high expression in non-human primates(NHP) and high tolerability in rats via intravenous (i.v.) administration.,

[Methods] LNPs encapsulating human erythropoietin (hEPO) mRNA were prepared using **FL-1245T** or benchmark ionizable lipid (Lipid 5). Plasma hEPO expression levels were evaluated in NHP following i.v. dosing. Plasma levels of liver enzymes and cytokines were evaluated in rats.

[Results] **FL-1245T**, our novel ionizable lipid, demonstrated higher hEPO expression than Lipid 5 in NHP. **FL-1245T** also demonstrated greater tolerability than our previously reported lipids and comparable tolerability to Lipid 5 at high doses (3mg/kg, 6mg/kg) administration in rats.

FUJIFILM-Gr mRNA-LNP Business

Formulation prototype Clinical trial manufacturing

Research Preclinical Ph 1 Ph 2 Ph 3 Commercial

• FUJIFILM ionizable lipids • Scale-up & manufacturing for toxic study • Commercial manufacturing

Lipid Licensing

Discovery Phase:

- Provide **FUJIFILM**'s proprietary ionizable lipids depending on client's purposes

Track Record : **FL-0445**

- ✓ Used as client's vaccine in clinical trial Ph3
- ✓ GMP manufacture is available

Kanagawa/JP Lab

CDMO Service

Research Phase:

- mRNA design and synthesis
- mRNA-LNP formulation optimization

Development Phase:

- process and analytical methods development
- mRNA-LNP manufacture in GMP

Toyama/JP Factory

FUJIFILM Lipid List

Application	Lipid Name	GMP Mfg.	IP Filled
Vaccine/ IM	FL-0445	✓	✓
	FL-2266	✓	✓
Liver, IV	FL-1245T		✓
	FL-1207T		✓
Spleen, IV	FL-1030T		✓
	FL-1625T		✓
T cell mRNA Delivery	FL-1779T		✓
T cell DNA Deliv., Liver-detargeting	FL-1779T		✓

FUJIFILM


Lipid 5 Lipid X FL-1030T

Fig. 1 Fig. 2 Fig. 3



Li: liver, Lu: lung, S: spleen K: kidney, H: heart

LNPs encapsulating fluc-mRNA were administered to ICR mice at a dose of 0.2mg/kg (i.v.) and fluc luminescence was imaged and quantified by IVIS. (Top: 5hr, Bottom: 6hr)

Result 1: NHP study revealed that FL-1245T achieved high protein expression



Result 2: Rat study revealed that FL-1245T is highly tolerable



Conclusions and Discussion

➤Our novel ionizable lipid **FL-1245T** demonstrated

- high protein expression following a low-dose i.v. infusion in monkeys.
- high tolerability following a high-dose i.v. administration to rats.

➤To expand the potential of mRNA-LNP therapeutics:

- understanding the mechanisms of liver toxicity caused by ionizable lipids would be beneficial for further improving the design of novel ionizable lipids from both efficacy and safety perspectives.

Please Visit Our Booth No. 8 !

& other posters No.35, 117

Our Website & LinkedIn:

<https://www.fujifilm.com/jp/ja/business/cdmo/medical/liposome>

CDMO website LinkedIn